

Orbi-Turn Patient Turner



- The Orbi-turn has been designed for carers to facilitate standing and seat-to-seat transfers
- It allows a patient with moderate upper and lower limb strength to participate in functional transfers whilst ensuring good working postures for carers and therapists
- Safe Working Load 190 kg

Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

ndis

Item Code:

05_123106154_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

| Code | Display Name | RRP | Box Quantity |
|-----------|-----------------------------------|------|--------------|
| TSE672935 | Orbi-Turn Standing Patient Turner | 1519 | 1 |