

## Oxford Elevate Standing Lifter



- The Elevate is a high performance active lift designed to improve the lifting experience for both the caregiver and the patient
- It is both compact and sturdy with a safe working load rating of 200 kgs
- Integrated digital weigh-scale assist with patient monitoring, designed with intuitive controls that require very little product training for caregivers
- The foot-tray is angled and utilises a lightweight plastic cover that can be easily removed for cleaning. In addition, by removing the foot-tray, the Elevate allows gait training for rehabilitation of patients
- The sculpted kneepad can be easily adjusted to improve patient comfort
- Item Weight: 62 kg



**ndis**

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 2 Years



**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Electric Leg	Safe Working Load	Material
LSS377302	Elevate Patient Lifter	1100 mm	630 mm	1130 - 1600 mm	945 - 1640 mm	510 - 870 mm	True	200 kg	Aluminium