


Aspire Low Friction Sling with Head Support - Parasilk



- Australian standards tested and registered with the TGA
- Parasilk material, with low friction properties, for maximum patient comfort during transfers and easy positioning
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- Durable - loop webbing has a breaking strain in excess of one tonne and is designed so that it will not tangle or “cork-screw” during cleaning
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- Parasilk leg strap lining for easy positioning
- For FIM score 1 or 2

 Contract approved

ndis

Item Code:

05_123603166_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Safe Working Load	Material
LSS390880	Deluxe Low Friction Sling - Parasilk - Small	499	1	200 kg	Parasilk
LSS390885	Deluxe Low Friction Sling - Parasilk - Medium	499	1	200 kg	Parasilk

Code	Display Name	RRP	Box Quantity	Safe Working Load	Material
LSS390890	Deluxe Low Friction Sling - Parasilk - Large	499	1	300 kg	Parasilk
LSS390895	Deluxe Low Friction Sling - Parasilk - Extra Large	499	1	300 kg	Parasilk