

## Feeding Cup



- A clear plastic beaker, designed for easy gripping, which has a vented lid with a drinking spout
- Available with 4 mm or 8 mm spout hole
- Sold in pairs
- Weight: 80 - 100 g

### Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

### ndis

#### Item Code:

05\_150900319\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Capacity
DLK276550	Feeding Cup - 4mm Spout	200 ml
DLK276560	Feeding Cup - 8mm Spout	200 ml