

Freedom Toilet Seat Raiser - HBA444



- Padded armrests can swing up and out of the way or be easily removed for easy cleaning
- Entire seat flips up so the bowl can be accessed without the seat (except on lowest height setting)
- Height can be adjusted to 5, 10 or 15 cm to ensure hips remain above knee height
- Attaches to existing toilet seat mounting points

- ✔ Contracted DVA product
- ✔ Contract approved

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year



Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth
BTT164100	Freedom Toilet Seat Raiser	380 mm	50, 100 or 150 mm	136 kg	3.8 kg	430 mm