



iCare Bed



- The IC333 Homecare Bed is unmatched for function and durability, making the hospital look out of your room
- This four-function bed features dual lift motors totaling 12,000 newtons lifting power, unlimited under-bed clearance, and eight heavy-duty casters
- Hand Piece – Backlit for night time use. Easy press buttons. Wired stretch cord. Hook attachment for hanging to side of bed. Wipeable
- Mattress Surface – Strong metal mesh under mattress allowing increased airflow
- Wheels – 8 x heavy duty swivel castors with one lockable castor on each corner
- Head lift angle 62° and knee break angle 42°
- End and side mattress retainer bars
- Ultra low design reducing risk of falls
- The safe working load includes any model of iCare mattresses
- Headboard and mattress not included.

 Contracted DVA product

 Contract approved

ndis

Item Code:

05_181210201_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

IEC 60601-2-52:2009

Warranty Information

- 10 year Limited Warranty

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB045952	IC333 Bed Base - Onyx - Long Single	2150 mm	960 mm	220 - 660 mm	200 kg	200 kg
BEB045958	IC333 Bed Base - Onyx - King Single	2150 mm	1120 mm	220 - 660 mm	200 kg	200 kg
BEB045963	IC333 Bed Base - Onyx - Queen	2150 mm	1580 mm	220 - 660 mm	300 kg	200 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
BEB045965	IC333 Bed Base - Onyx - Long Double	2150 mm	1400 mm	220 - 660 mm	300 kg	200 kg
BEB045966	IC333 Bed Base - Stone - Long Single	2150 mm	960 mm	220 - 660 mm	200 kg	200 kg
BEB045967	IC333 Bed Base - Stone - King Single	2150 mm	1120 mm	220 - 660 mm	200 kg	200 kg
BEB045968	IC333 Bed Base - Stone - Long Double	2150 mm	1400 mm	220 - 660 mm	300 kg	200 kg
BEB045969	IC333 Bed Base - Stone - Queen	2150 mm	1580 mm	220 - 660 mm	300 kg	200 kg
BEB045670	IC333 Bed Base - Black Vinyl - Long Single	2030 mm	900 mm	660mm - 220mm	200 kg	-
BEB045951	IC333 Bed Base - Vinyl Ebony - King Single	2030 mm	1070 mm	660mm - 220mm	200 kg	-
BEB045963-S	IC333 Bed Base - Onyx - Queen Split	2030 mm	1520 mm	660mm - 220mm	200 kg	-
BEB045969-S	IC333 Bed Base - Stone - Queen Split	2030 mm	740mm (per side)	660mm - 220mm	200 kg	-