

The Anaesthetic Cart



- Adjustable push handle, quick and easy operation to enable users to set the push handle at the correct ergonomic height
- Recessed ABS formed plastic work surface
- 12 x tilt storage containers
- Adjustable intravenous pole
- Slide out shelf fitted with runners for smooth, stable operation
- Waste bin
- Sharps holder
- Lightweight composite aluminium/plastic panels
- Corner buffers
- Raised drawer handles designed for grip and easy cleaning
- Diameter 125 mm anti-static total lock castors
- All drawers fitted with soft close runners and damper - prevents drawers from opening while in motion
- Epoxy powder coated frame



Warranty Information

- 5 Years Frame

Code	Display Name	Total Length	Total Width	Total Height
WAT901170	Anaesthetic Cart (Blue) - 5 Drawer + S/Steel Top + Guardrail + Tilt Containers (17)	545 mm	760 mm	1580 mm