

Freedom Bariatric Walking Frame - RBE804



- Heavy duty strong steel construction
- PVC handgrips
- Height adjustable extension legs
- Replaceable skid resistant rubber tips
- Easy-fold mechanism for storage and transportation
- Folded width: 95 mm

✔ Contracted DVA product

✔ Contract approved

Compliance

ISO 11199-2:2005

Warranty Information

- 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight
WAF752162	Freedom Bariatric Walking Frame	460 mm	740 mm	760 - 1010 mm	227 kg	Steel	3 kg