

## 4H175E Electric Paediatric Cot



This cot is thoughtfully designed for paediatric care and caregiver convenience, offering features that improve safety and efficiency in hospital settings.



- **Easy-to-Operate Side Rails:** The side rails are easy to operate with one hand, enhancing patient safety and allowing quick adjustments, even when multitasking.
- **All-Around Access:** The open base design allows access from all sides, making patient care easier and reducing the risk of injury.
- **Height Customisation:** The electric height adjustment makes it easy to position the cot at the right level for both patient comfort and caregiver convenience.
- **Emergency-Ready Positioning:** The cot can quickly tilt up to +12 degrees for emergency situations. For stability, use the manual release lever to return the cot to a flat position instantly.
- **Integrated Battery Backup System:** The battery backup ensures that height and tilt adjustments can be made during power outages, providing consistent care in all situations. It's not intended for long-term cordless use.

### Warranty Information

- 5 Years - Frame
- 2 Years - Electronics
- 1 Year - Hand Control & Cables

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Colour
WAT830560	4H175E Electric Paediatric Cot - White/Green	1410 mm	880 mm	1490 mm	50kg	White/Green
WAT830570	4H175E Electric Paediatric Cot - Custom Rainbow	1410 mm	880 mm	1490 mm	50kg	Custom Colour
BEM046980	Mattress - 100mm	-	-	-	-	-
WAT740000	Fixed Height IV Pole	-	-	-	-	-
WAT740002	Fixed Height IV Pole - w/25mm round tube to mount Pump	-	-	-	-	-
WAT740004	Height Adjustable IV Pole with Pump Mount	-	-	-	-	-