

JAY - Lite



- An ideal off-the-shelf complement to ultra lightweight wheelchairs
- Constructed with lightweight 3DX Microclimatic spacer fabric and Airflow base to provide excellent microclimatic heat and moisture dissipation
- The OptiWell pelvic loading area helps to achieve high levels of pressure distribution and stability in an extremely lightweight cushion
- Available with positioning inserts for those with postural needs (please ask our Aidacare C.A.S.E consultants for help selecting these)
- Available in sizes 14 x 14" to 24 x 24"
- All models include: Incontinence resistance cover

- ✔ Contracted DVA product
- ✔ Contract approved

ndis

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year Cushion
- 6 Months Cover



Compliance

No applicable compliance standard

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC571913	JAY Lite - 16 x 18in	2.25 in	113 kg	Medium Risk	Spacer Fabric
PTC585157	JAY Lite - 18 x 16in	2.25 in	113 kg	Medium Risk	Spacer Fabric