

## ISA XPLUS Stander



- Three foot plate heights available with removable heel block and round shape for toilet approach. Intermediate foot plate supplied with lifter
- The push bar and foot pedal have an ergonomic shape for easy lifter operation
- Sling hooks are flexible and secure for a quick, straightforward sling attachment
- Leg band padded for extra comfort. Easy to exchange and washable
- Curved legs have an ergonomic shape for easy manoeuvring around chairs and toilet bases
- Lockable rear casters are added security and safety for the patient and user
- Removable battery has a second battery charging to reduce lifter down time



Code	Display Name	Total Length	Total Width	Safe Working Load	Lift Range
LSS378050	Invacare Lifters - ISA XPLUS Stander	1280 mm	640 - 1020 mm	200 kg	505 - 675 mm