

Mortuary Spreader Bar 363kg



- To be used in conjunction with the Morgue Sling (4 bands) to allow optimal transfers in morgue applications
- The carry bar has a SWL of 363kg
- Accessory is suitable for use with Fixed Ceiling Lifts
- Compatible with Liko patient lifters & Liko ceiling hoists, Prism ceiling hoists, Aidacare A320, A205 and 200 patient lifters



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

sova / dignity in motion

| Code | Display Name | Safe Working Load |
|-----------|--|-------------------|
| LIC354197 | Mortuary Spreader Bar Lifting strap - Short - 2 eyelets - HDPE Plastic - 50mm wide - | |
| LIC362300 | Mortuary Spreader Bar - 363 kg | 363 kg |
| LIC354198 | Mortuary Spreader Bar Lifting strap - Short - 4 eyelets - HDPE Plastic - 50mm wide - | |