

Tilley Wingback Chair



- A practical chair featuring a contemporary take on a classic wingback
- The timeless sophistication of the straight, clean lines will complement a variety of interior styles
- The generous seating area looks and feels inviting
- Tapered legs in a seasoned oak finish provide a subtle finish
- Foam cushioning and quality construction ensures this armchair will withstand daily use
- A thoughtful addition to a bedroom or lounge area
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052240	Tilley Wingback Chair	750mm	700mm	975mm