

## Liko MasterVest (Mod. 64)



- Versatile training aid that facilitates toilet visits, standing training, gait training, mobility training and placing patients into standing positions
- Reduces strain on caregivers while helping them deliver physical training
- Transfers to and from the toilet can be accomplished easily and smoothly, reducing strain
- Lifts around the upper body and partly under the arms, and is best suited for patients that can support some of their body weight
- Optional leg harness for improved support [Liko Mastervest Walking Sling Leg Harness](#)
- Available in Polyester

### Warranty Information

- 1 Year



### ndis

#### Item Code:

05\_123603166\_0105\_1\_2

#### Hoist –

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394260	MasterVest Walking Sling - Front Close Mod 64 - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394270	MasterVest Walking Sling - Front Close Mod 64 - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394280	MasterVest Walking Sling - Front Close Mod 64 - Polyester - LARGE	200 kg	Polyester (Loops)
LSS394290	MasterVest Walking Sling - Front Close Mod 64 - Polyester - XLARGE	200 kg	Polyester (Loops)
LSS394300	MasterVest Walking Sling - Front Close Mod 64 - Polyester - XXLARGE	200 kg	Polyester (Loops)