

Oxford Standing Sling



- Attaches to hoist to enable a carer to lift and stand a less mobile user
- High density foam maximises user comfort and support
- Quick and easy to fit adjustable velcro waist band
- Colour coded attachment and positioning straps
- Designed for use with standing hoists
- FIM Scores: 3 and 4

Contracted DVA product

Contract approved

ndis

Item Code:

05_123603162_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS371951	Standing Sling - Medium	200 kg	Polyester
LSS371950	Standing Sling - Large	200 kg	Polyester