

## Panamatic Optima 2



- The Panamatic Optima 2 - Flusher Thermal Sanitiser (the hot & cold water cousin of the midi) has raised the standard to help you meet the ongoing hygiene challenges in the Hospital and Aged Care environments
- Improved Wash Performance (Hot & cold water) - An increased number of more efficient cleaning nozzles and optimum water jet pressure result in performance levels that make the Optima 2 outstanding for its size
- Reduced Cycle-Time - Better washing and temperature control has enabled improved overall cycle-time for quicker turnaround
- More Efficient Water and Power - With improved performance and reduced cycle time, both power and water consumption have been reduced. Cold water flush / wash - steam cycle - cool-down rinse
- Hands-free Opening and Cycle Start - A foot operated opening mechanism and hands-free close and cycle start means direct contact with the machine has been reduced to a minimum, helping to reduce potential cross infection
- MicrobeSafe Antimicrobial Lid - DDC™ unique lid incorporates Silver Nanotechnology to help stop bacterial growth and reduce cross infection
- Large Chamber Capacity and Easier Loading - The wash chamber and larger opening allow easier loading of a wide range of utensils. The chamber can process 2 bedpans and 2 urine bottles per load or various combinations
- Electrical supply: Single Phase 230 V, 13 AMP
- Electrical output/input: Motor - 0.75 KW, Heater - 2.5KW
- Water supply: 15 mm Hot & Cold 1 - 6 bar
- Soil connection: 110 mm connections supplied
- Operation cycle time: Approx 5 mins / 30 ml EcoWash
- Steam disinfection: Exceeds Ao 700, 85° x 75 secs



## Warranty Information

- 1 Year

Code	Display Name
UPM100210	DDC Bedpan Washer Sanitiser - Optima 2