


Gel Armrest Pads



- Ultimate gel comfort and protection for the elbows and wrists
- Distributes pressure efficiently to help avoid harmful pressure-related problems
- Durable neoprene nylon surface for easy cleaning

 Contract approved

Compliance

No applicable compliance standard

Warranty Information

- 1 Year

| Code | Display Name | Total Length | Total Width |
|-----------|-----------------------------|--------------|-------------|
| PTA508350 | Gel Armrest Pads 3.5 x 10in | 254 mm | 88.9 mm |
| PTA508400 | Gel Armrest Pads 3.5 x 12in | 304.8 mm | 88.9 mm |
| PTA508425 | Gel Armrest Pads 3.5 x 14in | 355.6 mm | 88.9 mm |
| PTA508450 | Gel Armrest Pads 3.5 x 16in | 406.4 mm | 88.9 mm |