

Aspire Bariatric Shower Commode



- Super Heavy Duty Bariatric design with 790 mm seat and 350 kg SWL
- Anti-bacterial frame coating for improved infection control
- Removable arms to assist with positioning and allow lateral slide transfers
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- Padded seat and arms for maximum comfort
- 4 x lockable castors maximise safety by providing carer access to a brake at all angles around the commode
- Comes standard with pan and carrier
- Overall Height: 990 mm
- Overall Width: 900 mm
- Backrest Height: 430 mm
- Armrest Height: 170 mm
- Anterior Seat Height: 570 mm
- Posterior Seat Height: 570 mm

Prescribing Therapists: Click here to download our [Clinical Assistance Guide](#)

Warranty Information

- Corrosion: 5 Years
- Components: 1 Year

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_091203054_0103_1_2

Mobile Shower Commode –

Composite:05_413_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

Accessories

Click on the cart icon to display the range of optional accessories

Code	Display Name	Seat Width	Safe Working Load
BTC061200	Bariatric Shower Commode	790 mm	350 kg