

## Aspire Hybrid Patient Sling



- Hybrid sling facilitates both Hygiene and General Purpose manual handling
- Intuitive shape allows ease of access for toileting and dressing
- Full coverage shoulder and head support
- Security belt and padded leg straps support the trunk and thighs
- Durable loop webbing designed not to tangle or 'cork-screw' during cleaning
- Embroidered unique serial numbers for traceability
- Service record labels to easily record and identify service history
- Machine washable up to 85°



### Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Safe Working Load
LSS390830	Aspire Hybrid Patient Sling - Small - 200kg	475	1	200 kg
LSS390835	Aspire Hybrid Patient Sling - Medium - 200kg	475	1	200 kg
LSS390840	Aspire Hybrid Patient Sling - Large - 300kg	475	1	300 kg
LSS390845	Aspire Hybrid Patient Sling - Extra Large - 300kg	475	1	300 kg