

Trekker Gait Trainer



- Cardiopulmonary strengthening and conditioning
- Reduces flexion contractures to enhance movement
- Increases spatial awareness and head control to achieve independent ambulation
- Teaches individual to take steps without assistance
- Can be used in the anterior and posterior position providing flexibility for the user
- Directional lever allows both forward and reverse mobility or only forward mobility
- Available in 2 sizes - Specify size on request
- Accessories available on request



ndis

Item Code:

05_120609099_0105_1_2

Walking Aids -

Composite:05_416_0105_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Seat Depth	Seat Width	Safe Working Load
PED480820	Trekker Gait Trainer	790 - 940 mm	620 - 670 mm	68 - 90 kg