

Aspire Classic Bedside Commode



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Raised Padded Armrests
- Easy Clean Vinyl Surfaces
- Height Adjustable
- Removable, Padded Seat
- Non-Slip Rubber Tips
- Toilet Seat to Floor Height 475 - 600 mm

Contracted DVA product

Contract approved

ndis

Item Code:

05_091203048_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BEC045005	Aspire Classic Bedside Commode - Champagne	440 mm	450 mm	500 - 625 mm	125 kg	5.5 kg
BEC045006	Aspire Classic Bedside Commode - Slate	440 mm	450 mm	500 - 625 mm	125 kg	5.5 kg