

Start JNR



- Lightweight paediatric chair
- Full range of accessories to match the needs of the child
- Now available in a range of colours



ndis

Item Code:

05_122203114_0105_1_2

Wheelchair Manual –

Composite:05_417_0105_1_2

Click [here](#) for more information



Code	Display Name	Seat Depth	Seat Width	Total Width	Backrest Height	Seat To Floor Height	Safe Working Load	Unit Weight
MWC427410	Start Junior 30 cm (as per script)	24 - 40 cm	28 - 38 cm	-	-	395 - 520 mm	90 kg	13.5 kg
MWC427420	Start Junior 35 cm (as per script)	-	-	355mm	310 - 370 mm	-	90 kg	-