

Oxford Ascend Standing Lifter



- With it's great range of movement, the Ascend can transfer a patient from a low-seated position to fully extended with ease
- The Ascend is compact and easy to manoeuvre in confined spaces such as bathrooms and has many key design features that make it a truly flexible piece of equipment
- Its active lifting motion has been designed to encourage user participation, which ultimately promotes patient independence and well-being
- The kneepad can be adjusted within a 150 mm operating range to improve patient comfort during repositioning
- The removable foot tray allows the Ascend to be used as a walking aid - great for patients that require gait-training
- Ergonomic foot pad to assist with initiating movement
- Careful attention to product design ensures correct ergonomic usage, protecting both the patient and carer from injury at all times
- Item Weight: 47 kg



ndis

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 2 Years

| Code | Display Name | Total Length | Total Width | Total Height | Lift Range | Leg Gait Range | Electric Leg | Safe Working Load | Material |
|-----------|-----------------------|--------------|-------------|----------------|---------------|----------------|--------------|-------------------|-----------|
| LSS377350 | Ascend Patient Lifter | 1100 mm | 630 mm | 1070 - 1710 mm | 990 - 1695 mm | 510 - 885 mm | True | 170 kg | Aluminium |