

Hensinger Head Support



- Hensinger features a soft foam wrap around the collar to support the head by using base support areas at the occiput and jaw line
- Hensinger head support must be used in conjunction with a chest support
- This head support promotes proper head positioning resulting in improved breathing, swallowing and simplified feeding
- Hensinger enhances more social interaction and awareness during daily activities
- Revolutionary in the world of cervical collars, this soft foam wraparound collar is designed to brace the head by the jawline
- The Hensinger Head support can help maintain mid-line head positioning. Proper head positioning facilitates improved breathing, swallowing, and simplified feeding
- It also aids in social interaction and awareness during classroom activities
- The Hensinger is built with a strong yet pliable reinforcement material around the outside of the support
- A simple Quick-Release buckle allows easy on/off
- The Hensinger Head Support can be used with a chest support or mounted to a seating system
- Soft covers are included
- **Sizing guide:**
 - Length measurement = neck circumference
 - Height is shoulder to bottom of ear, allowing 13 mm between earlobe and top of head

✓ Contract approved

Support @ Home

Function: Body-worn assistive products for body protection

Category: Under advice

Click [here](#) for more information

Code	Display Name	Total Length	Total Width
FBS279900	Hensinger Head Support - X-Small	228 mm	25 mm
FBS279901	Hensinger Head Support - Small	304 mm	38 mm
FBS279902	Hensinger Head Support - Medium	355 mm	50 mm
FBS279903	Hensinger Head Support - Large	406 mm	63 mm