

## Heelift Ultra Suspension



- Ideal for patients with foot-related mobility issues
- Prevents pressure ulcer and provides comfort to users with existing pressure ulcers
- Smooth, durable exterior fabric prevents boot from catching on sheets and blankets
- Forefoot strap keeps boot in place and provides support
- Adjustable, soft rounded fastening straps
- Open cell construction provide additional air circulation
- Friction-free tricot backing increases patient mobility
- Deluxe exterior stitching provides durable boot structure
- Easy-to-clean, lightweight fabric with inner lining that holds the leg in place
- Machine washable



Code	Display Name
PTA509085	Heel Protector - Heelift Ultra Suspension - Standard
PTA509090	Heel Protector - Heelift Ultra Suspension - Petite
PTA509095	Heel Protector - Heelift Ultra Suspension - Bariatric