

## MOTomed Loop.LA Pro



**RECK**  
MOTomed®

The MOTomed Professional Series are tailor-made devices for versatile use in clinics and rehabilitation facilities. The basic equipment of MOTomed loop.la prof includes:

- The LA comes with both upper limb and lower limb trainers installed
- Plastic-coated safety foot shells
- Leg guides with plastic-coated calf shells and disinfectable fastening straps
- QuickFix foot fastening system
- Plastic-coated handles with quick release system
- Gas spring support for height adjustment
- Pedal radius quick adjustment (leg trainer, pair)
- Special stand with reduced stand height
- Activation of serial interface
- Serial port activated: trainer can be upgraded to an FES Cycle with the HASOMED RehaMove2 'Functional Electrical Stimulator'
- Hand fixation by wrist cuff (piece)



### Accessories

Click on the cart icon to display the range of optional accessories

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101120	MOTomed Loop.LA Pro - Serial Interface, Special Stand, Hand Fix	70 cm	60 cm	107 - 122 cm	135 kg	39 kg