

TENA Slip



- Specially designed for those care recipients with greater needs and who could be either bed-ridden or over-active
- All-in-One product and features contoured, elasticized edges that ensure maximum security and comfort



ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

No applicable compliance standard



Code	Display Name	Capacity
CON242000	TENA Slip - Extra Small - 30 Pack	1410 ml
CON242300	TENA Slip - Maxi Large - 9 Pack	2300 ml
CON242200	TENA Slip - Maxi Medium - 9 Pack	2000 ml
CON242100	TENA Slip - Maxi Small - 24 Pack	1200 ml
CON242600	TENA Slip - Plus Large - 12 Pack	1500 ml
CON242500	TENA Slip - Plus Medium - 12 Pack	1300 ml
CON242900	TENA Slip - Super Large - 10 Pack	1800 ml

Code	Display Name	Capacity
CON242800	TENA Slip - Super Medium - 10 Pack	1600 ml
CON242700	TENA Slip - Super Small - 28 Pack	950 ml