

Aidacare Aspire 200S Steel Standing Lifter



The Aidacare Aspire 200S is a robust and durable steel standing lifter capable in the most demanding facility environments

- **Bullhorn Lifting** Handles for patient stability and 6 sling attachment points
- **Ergonomic Handles** Provide multiple grip points to suit carers of different heights and strengths
- **Ergonomic Handset** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Powered Leg Function** Opens and closes the legs for lifting stability or accessibility
- **Sturdy steel frame** with industry leading Linak electronics and warranty provision
- **Safety** Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- **Adjustable Padded Legrests** Provide comfort and support during transfers & can be adjusted to suit the patient
- **Intelligent Diagnostics** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- **Tente Casters:** 2 x 75 mm front castors for maximum under bed clearance and accessibility. 2 x 75 mm rear locking castors for safety and manoeuvrability
- Item Weight: 44.6 kg
- Tested and compliant to AS/NZS ISO 10535:2011 (the Australian standard for hoists)

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_123604167_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- Frame: 5 Years
- Electronics: 3 Years
- Battery and Handset: 1 Year

| Code | Display Name | Total Length | Total Width | Total Height | Safe Working Load | Material | Unit Weight | Lift Range | Leg Gait Range | Electric Leg |
|-----------|--|--------------|---------------|--------------|-------------------|----------|-------------|---------------|----------------|--------------|
| LSS390650 | Aidacare Aspire 200S Standing Patient Lifter | 1150 mm | 656 - 1016 mm | 1790 mm | 200 kg | Steel | 44.6 kg | 850 - 1740 mm | 575 - 936 mm | True |