

Baby Bassinet Mobile



- Moulded acrylic bassinet
- 2 x moulded plastic utility trays
- Soft vinylclave covered mattress
- Bassinet can be elevated at either end

Warranty Information

- 3 Years



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
WAT829000	Baby Bassinet Mobile - Single Acrylic	780 mm	400 mm	990 mm
WAT829600	Baby Bassinet Mobile - Single S/Steel	780 mm	400 mm	990 mm