

Medical Surgical Masks - Disposable (Level 2)



- **Level 2:** > 98% Bacterial Filtration Efficacy.
- Standards Compliance: **EN14683 and ISO22609**
- CE Compliant
- Multi-layered Non-Woven Construction (Spunbond)
- High Frequency Welding for strongest bonding and elimination of holes and defects
- Bacterial Filtration Efficiency: 99.7%
- Particle Filtration Efficiency: 98.5%
- Particle size: 0.1 μm

These masks are disposable and intended as single use (up to 1 day if not soiled)

Box & Order Quantity

Packaged: 50pc/bag; 40bags/box

MOQ - Individual Orders = 1 Box = 2,000 units



General Order Information

Lead Time 5 days + 5-days

Shipping = 10days

Sea Freight = 3 weeks

MOQ TOTAL ORDER - 100,000pc

Payment terms are 100% prepayment