


Aspire Autoclavable Wash Bowls



Our range of Wash Bowls are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification
- *Length measurement refers to the diameter

 Contract approved

Compliance

No applicable compliance standard



Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Length	Capacity
WAT769800	Wash Bowl - Small	-	-	140 mm*	600 ml
WAT770000	Wash Bowl - Medium	-	-	210 mm*	1500 ml
WAT761200	Wash Bowl - Large	28	1	305 mm*	2400 ml
WAT769600	Wash Bowl - Extra Large	36	1	345 mm*	5500 ml