

## MOTomed Wheelchair Stabiliser



**RECK**  
MOTomed®

- Recommended for severe spasticity
- Supports the wheelchair: prevents tipping and slipping backwards
- For all common wheelchairs with push handle
- With padded support for the wheelchair push handles
- Width of the handle support: 580 mm Children's stabiliser especially designed for children's wheelchairs with push handle: Width of the handle support: 330 mm



Code	Display Name
FPM102170	MOTomed Wheelchair Stabiliser
FPM102175	MOTomed Wheelchair Stabiliser - Children