


## K Care URETHANE Commode Seat Padded - Open



- High quality urethane padded seat with an open front for personal cleaning and toileting access
- Available in blue
- Supplied with plastic seat runners for easy attachment and removal during cleaning of the seat and commode
- Seat Height: 40 mm

 Contract approved

**Compliance**  
AS/NZS 3973:2009

### Warranty Information

- 1 Year

| Code      | Display Name                          | Seat Depth | Seat Width |
|-----------|---------------------------------------|------------|------------|
| BTC075700 | Open Front Commode Seat - 445 mm Wide | 430 mm     | 445 mm     |