

Aidacare FL250 Floorline Bed - Telescopic Rails



- Adjustable height rails to assist in the prevention of falls as well as manoeuvring.
- Includes both left/right side rails, and bed height limiting handset.
- Various configurations include 1/3, 2/3 individually or together for full length
- *Bed ends sold separately
- Please note: The minimum bed height changes to 275mm when Telescopic Rails are installed



| Code | Display Name | Colour |
|-----------|-----------------------------------------------------------------------------|---------------|
| BEB048800 | FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Rail - Sublime Teak | Sublime Teak |
| BEB048805 | FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Rail - Milano Walnut | Milano Walnut |
| BEB048810 | FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Decorative Trim - Sublime Teak | Sublime Teak |
| BEB048815 | FL250 Telescopic Rails - 2x1/3 Rail + 2x2/3 Decorative Trim - Milano Walnut | Milano Walnut |
| BEB048820 | FL250 Telescopic Rails - 2x2/3 Rail + 2x1/3 Decorative Trim - Sublime Teak | Sublime Teak |
| BEB048825 | FL250 Telescopic Rails - 2x2/3 Rail + 2x1/3 Decorative Trim - Milano Walnut | Milano Walnut |
| BEB048830 | Bed End Set - FL250 - Single - ComoR (Sublime Teak) - For Telescopic Rails | Sublime Teak |

| Code | Display Name | Colour |
|-----------|-----------------------------------------------------------------------------------|---------------|
| BEB048835 | Bed End Set - FL250 - Single - ComoR (Milano Walnut) - For Telescopic Rails | Milano Walnut |
| BEB048840 | Bed End Set - FL250 - King Single - ComoR (Sublime Teak) - For Telescopic Rails | Sublime Teak |
| BEB048845 | Bed End Set - FL250 - King Single - Como R (Milano Walnut) - For Telescopic Rails | Milano Walnut |