

## MacIntosh Waterproof Protector



- Aimed at increasing the life span and decreasing the drying time of the bed pad
- Comes standard with tuckins to secure firmly on the bed
- Ideal for use with CON225008, CON228600, CON230800 & CON231100



**ndis**

**Item Code:**

03\_092489060\_0103\_1\_1

Click [here](#) for more information

### Warranty Information

- 1 Year

### Compliance

No applicable compliance standard



Code	Display Name	RRP	Box Quantity	Total Length	Total Width
CON229100	MacIntosh Waterproof Protector	52	1	900 mm	1000 mm