

Thera-Band Soft Weights

- Provide more versatility for strength training and rehabilitation exercises
- Hand-held, ball-shaped isotonic weights



Code	Display Name	Unit Weight	Colour
FPF290381	Thera-Band Soft Weights - Tan 0.5kg	0.5 kg	Tan
FPF290382	Thera-Band Soft Weights - Yellow 1.0kg	1.0 kg	Yellow
FPF290383	Thera-Band Soft Weights - Red 1.5kg	1.5 kg	Red
FPF290384	Thera-Band Soft Weights - Green 2.0kg	2.0 kg	Green
FPF290385	Thera-Band Soft Weights - Blue 2.5kg	2.5 kg	Blue
FPF290386	Thera-Band Soft Weights - Black 3.0kg	3.0 kg	Black