


## Liko Sabina II EE Lifter



- Conventional patient lifting to a limited extent, lifting in the sitting position
- Powerful lifting motor ensures the smooth raising and lowering of patients
- Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility
- The lower leg support and the foot rest can easily be removed for cleaning
- Adjustable height settings accommodate a variety of different patient sizes
- Sling bars are sold separately to give customer the choice of best one for their needs. We recommend either the Sabina II 350 SlingBar (LSS393270) or Sabina II Comfort Sling Bar (LSS393280)

### Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

 Contract approved

**ndis**

**Item Code:**

05\_123604167\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393250	Liko Sabina II EE Lifter	1060 mm	695 - 1115 mm	900 - 1000 mm	200 kg	Steel