

Contego Floorline Folding Bed



- Floor bed lowers to within 110 mm off the floor for safety and access in aged care
- Fully height adjustable with Trendelenburg tilt (15°) for drainage and positioning
- 75° Backrest and 30° Kneebend
- Lightweight, dismantle frame for easy transportation
- Integrated cable management to suit all alternating air mattress systems
- Removable head and footboards
- 100 mm total castors
- Standard single bed or king single bed extension available
- Mattress retainers are included
- Recommended mattresses: Lifecomfort Deluxe regular (BEM046500) or Lifecomfort Deluxe King Single (BEM046550)



Compliance

IEC 60601-2-52:2009

Warranty Information

- 10 Years Frame
- 3 Years Electrics