

Queens Cutlery



- A wide range of stainless steel utensils mounted into a built-up plastic handle to help those with weak or restricted grip, the handle is cylindrical to assist with gripping
- The items are hygienically sealed, easy to clean and are dishwasher safe
- Handle Length: 98 mm
- Handle Diameter: 30 mm
- Product Weight: 45 g

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name
DLK275760	Knife
DLK275770	Fork
DLK275780	Dessert Spoon
DLK275790	Junior Spoon
DLK275870	Cutlery Set