

McDavid A310 Athletic Support



- A multi-purpose supporter featuring a power-knit pouch with Intera for rapid drying and cooling
- Provides greater resistance to bacteria and mildew and the odours and skin irritation they produce
- More comfortable than cotton or synthetics

- ✔ Contracted DVA product
- ✔ Contract approved

Code	Display Name	Waist/Hip Circumference
FBS3016184	McDavid - Scrotal Support - Athletic Belt - XL	1120 - 1280 mm
FBS3016183	McDavid - Scrotal Support - Athletic Belt - Lge	970 - 1120 mm
FBS301618	McDavid - Scrotal Support - Athletic Belt - Sml	660 - 820 mm
FBS3016182	McDavid - Scrotal Support - Athletic Belt - Med	820 - 970 mm