

MOTOmed Leg Guides With Calf Shells



RECK
MOTOmed®

- For secure support and guidance of the legs
- Spring-loaded calf shells nestle elastically to the body
- Comfortable fixation by means of fastening straps
- Articulated and height-adjustable (9.0 cm to 16.0 cm)
- Removable calf shells, fold-down guide rail (lower part)
- Only available in combination with Soft Grip safety foot shells (order no. 265.400)



Code	Display Name
FPM102280	MOTOmed Leg Guides With Calf Shells (pair)