


Aidacare Aspire Inline Weigh Scale



- Designed to be fitted on Aidacare Aspire Patient Lifters
- Allows for an accurate weight measurement to be taken by a carer when a patient is suspended in a sling
- Large weight capacity of up to 320kg
- Easy to read LCD display

 Contract approved

Compliance

AS/NZS ISO 10535:2011



| Code | Display Name | Safe Working Load |
|-----------|----------------------------|-------------------|
| LSS390680 | Aspire In-Line Weigh Scale | 320 kg |