


Aspire Active Air Acute 8 Mattress Replacement



- The Active Air Acute 8 is specifically designed to meet the demands of Critical Care environments, maximising patient comfort and therapeutic outcomes
- 1 in 3 Cell Alternation Cycle with 24 Narrow Profile cells – Improves patient comfort by reducing space between cells and increasing patient support surface area
- Focused Heel Care – Removable heel cells facilitate rapid response to changes in heel therapy needs
- Dual Compressor Pump – Reduces noise and efficiently manages patient therapy at highest patient weight loads (250 kg)
- Hospital Grade Cover – Thicker weight polyurethane fabric with increased bleach resistance to handle modern hospital infection control standards. Suitable for Sodium Hydrochloride solutions
- Lateral turning through integrated air cells, allowing the mattress to turn a patient without the need to add an after-market turning device
- Drastically reduces the moving and handling stresses on Care Staff by rotating the patient up to 40 degrees. This technology is ideal for:
 - Alleviating asymmetric or laterally located pressure areas
 - Simplified in-bed bathing
 - Rapid changing of linen
 - Perfect positioning of manual handling slings and slide sheets

 Contract approved

Compliance

IEC 60601-1-2

Warranty Information

- Cover: 2 Years
- Air Mattress: 2 Years
- Pump: 2 Years

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Safe Working Load	Pressure Care Rating
PTM639400T	Active Air ACUTE 8 w/ turnassist – Single	4799	1	2000 mm	900 mm	200 mm	30 kg - 250 kg	Very High Risk
PTM639410T	Active Air ACUTE 8 w/ turnassist – King Single	6169	1	2000 mm	1050 mm	200 mm	30 kg - 250 kg	Very High Risk