


## Aspire Deluxe Bedside Commode



- Allows for toileting in the bedroom for users unable to travel between the bedroom and bathroom
- Available in Champagne or Slate Vinyl
- Easy Clean Vinyl Surfaces
- Removable Bowl with Lid & Carry Handle
- Height Adjustable
- Raised Padded Armrests
- Fold Away, Padded Seat
- Non-Slip Tips
- Toilet Seat to Floor Height 435 - 585 mm

 Contracted DVA product

 Contract approved

### ndis

**Item Code:**

05\_091203048\_0103\_1\_2

**Toilet Aids – Composite:**

05\_414\_0103\_1\_2

Click [here](#) for more information

## Warranty Information

- 1 Year



### Compliance

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Floor Height	Safe Working Load	Unit Weight
BEC045115	Aspire Deluxe Bedside Commode - Champagne	450 mm	460 mm	480 - 630 mm	160 kg	12 kg
BEC045116	Aspire Deluxe Bedside Commode - Slate	450 mm	460 mm	480 - 630 mm	160 kg	12 kg