

Smart Check Device By ROHO



permobil

ROHO

- Smart Check® by ROHO® is a cushion feedback system that is intended to assist the user by saving an inflation range during cushion setup. During subsequent inflation level checks, Smart Check is intended to indicate whether adjustments are needed to remain in the saved range
- Smart Check must be connected to a ROHO® Single Compartment Cushion with Sensor Ready® Technology
- Works exclusively with ROHO Sensor Ready HIGH PROFILE® and MID PROFILE Single Compartment Cushions
- Press-to-check inflation monitor
- Stores individual's inflation range
- Unlimited frequency to check / ability to check whenever & wherever
- Affordable, portable and durable
- Two setup modes: Smart Setup and Optional Setup
- Ability to reset – can store new inflation level for individual's physiological changes over time
- No weight limit – when Sensor Ready Cushion is properly sized to the individual
- Easy to read faceplate with ergonomic shape
- Easy to connect and disconnect - detachable for storage
- Lightweight and reliable with automatic power off
- Batteries available globally at retail outlets
- Long battery life with low battery indicator
- Multiple cover options for Sensor Ready Cushions
- Easy to clean and disinfect



Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name
MSP500488	Smart Check Device By ROHO