

## MOTomed Loop Light.LA



**RECK**  
MOTomed®

The MOTomed Loop brings movement therapy into the future. Through the new design and its digital intelligence the user can experience passive or active training in a completely new way. With the new and innovative height adjustment, users can adjust the position of the MOTomed loop without using tools. The Light model provides high quality primary care with all the common functions required from a movement therapy trainer - including Spasm Detection program for maximum user safety.



- The LA comes with both upper limb and lower limb trainers installed
- Soft Grip safety foot shells with fastening straps (pair, incl. pedal cranks)
- Soft Grip Handles (pair, incl. pedal cranks)
- Gas spring support for height adjustment
- Special stand with reduced stand height
- Leg guides with plastic-coated calf shells (pair)
- QuickFix foot fastening system (pair)
- Hand fixation by wrist cuff (piece)

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101140	MOTomed Loop.LA Light - Serial Interface, Gas Spring, Stand, Leg guides, Foot & Hand Fix	70 cm	60 cm	107 - 122 cm	135 kg	36 kg
FPM101141	MOTomed Loop LA Light Basic - Leg Guides	-	-	-	-	-