


Bed Pad - Smart Barrier



- No PVC Backing - only breathable, waterproof PU (polyurethane), 2 way stretch fabric as used on premium hospital mattresses
- 100% brushed polyester facing ideal for client protection and comfort

Warranty Information

- 1 Year

 Contract approved

ndis

Item Code:

03_092489060_0103_1_1

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Capacity
CON225002	Bed Pad - Smart Barrier - No Tuck Ins	900 mm	850 mm	2500 ml
CON225004	Bed Pad - Smart Barrier - With Tuck Ins	900 mm	850 mm	2500 ml