

Aspire ComfiMotion Care Bed



- Night Assist: Luminous underbed halo for safe entry and exit throughout the night
- Press and hold safety handset: Intuitive, large button interface enhances usability
- Enclosed frame design: Ensures optimal positioning when getting into bed and safely covers the moving base platform
- Safe-Lowering Technology: Ensures that back, leg and Hi-Lo functions only lower under the force of gravity reducing potential for accidents or injuries
- Box-on-box design facilitates easier entry and exit and allows true flush positioning of beds in companion or dual adjustable bed configurations.
- German engineered motors for smooth, quiet adjustment
- *Total Height with:
 - Low Profile Fixed Legs 280 – 580 mm (available upon request)
 - 75mm Standard Castors 350 - 650 mm
 - 125mm Accessibility Castors 400 - 700 mm (optional)Mattress choice may increase or decrease overall bed height
- *Mattress not included



Contracted DVA product



Contract approved

Compliance

IEC 60601-2-52:2009

Prescribing Therapists: Click here to download our
[Clinical Assistance Guide](#)

Warranty Information

- 10 Years Limited Warranty

| Code | Display Name | Total Length | Total Width | Total Height | Safe Working Load |
|-----------|---|--------------|-------------|--------------|-------------------|
| BEB046720 | Aspire ComfiMotion Care Bed - Portobello Grey - Long Single (Headboard Sep) | 2030 mm | 920 mm | 280 - 580 mm | 180 kg |
| BEB046730 | Aspire ComfiMotion Care Bed - Portobello Grey - King Single (Headboard Sep) | 2030 mm | 1070 mm | 280 - 580 mm | 180 kg |