


K Care Pressure Reducing Commode Seat



- Special foam design, engineered to improve the spread of pressure across the seat
- Compliments the user's pressure management routine
- Waterproof but breathable antimicrobial four-way stretch Dartex material to optimises user comfort
- Independently tested and proven to be effective in reducing pressure levels
- Seat height: 65 mm



Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

| Code | Display Name | Seat Depth | Seat Width |
|-----------|---|------------|------------|
| BTC082100 | Closed Front Commode Seat - Pressure Reducing - 440 mm Wide | 430 mm | 440 mm |
| BTC082000 | Open Front Commode Seat - Pressure Reducing - 440 mm Wide | 430 mm | 440 mm |