


## ETAC Cloo Fixed Height Toilet Seat Raiser with Arms



- The Cloo from Etac is a height adjustable toilet seat raiser with tiltable seat and detachable armrests for added support
- The arm supports fold up and away, or can be fixed in a folded down position for extra security if needed
- Both armrests are detachable, Cloo can also be used with only one arm support if space around the toilet is limited
- The seat can be easily adjusted to three different heights: low (60 mm), medium (100 mm) or high (140 mm) for maximum versatility
- Simple, elegant toilet seat raiser
- Detachable armrests either side of seat
- Mount to your existing toilet three height settings
- Easy to adjust height while mounted
- Can be set with forward facing tilt
- Detachable seat for easy cleaning
- Tool-less mounting to any toilet bowl

 Contract approved

Code	Display Name	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Width
BTT162030	ETAC Cloo Fixed Height Toilet Seat Raiser with Arms	600 mm	60 - 140 mm	150 kg	4.8 kg	400 mm